



Molly® 2.25

**Molly**®  
2.25 mL  
The ideal safety solution

Device Type	Disposable Auto Injector
Primary Container	2.25mL Long Pre-filled Syringe   RNS Compatible
Dosage Type	Fixed   Full Contents
Route of Administration	Subcutaneous
Needle Insertion	Manual
Activation Method	Needle Cover
Injection Feedback	Audible click at beginning & end of injection, viewing window

## Handling Instructions



Remove cap



Press & hold against injection site to activate



Discard in sharps container

## Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Manual needle insertion
- ✓ Audible, visual & tactile feedback
- ✓ One-handed operation
- ✓ Larger volume
- ✓ Permanently hidden needle
- ✓ Ergonomic pull- or twist-off cap
- ✓ Robust anti-roll design